

Temperature Control Effects in Packed Emitter Nanobore LC-MS/MS for Protein Digest Analysis

Mike S. Lee¹, Lee Sawdey², Amanda Berg², Gary A. Valaskovic²

¹Milestone Development Services, Newtown, PA ²New Objective, Inc., Woburn, MA

Overview

- Controlled heating of nanobore LC packed-tip emitter
- Demonstrated enhancement of chromatography at elevated temperatures
- Demonstrated operation of sub-2 μm material below 5,000 psi

Introduction

Column temperature control is a common parameter for method development in conventional (mm diameter) HPLC. Although there are numerous known advantages to elevated-temperature nanobore LC (improved chromatographic speed, peptide recovery) the widespread adoption of this platform has been limited by interface issues. The common packed-tip emitter used in proteomics applications offers challenges of integration of temperature control and application of high voltage. Commonly employed fused-silica tubing, however, has excellent thermal-transfer characteristics due to the high surface area/volume ratio. The effects of higher than ambient temperature (50-80 °C) on the separation and column back-pressure characteristics of a packed-tip emitter are investigated.

Materials & Methods

Instrumentation

- 3-D ion-trap mass spectrometer (LCQ Deca™; Thermo Fisher)
- Customized nanospray source equipped with ceramic cartridge heater assembly (Digital PicoView™)
- Direct flow capillary-LC pump (Agilent®)
- nano LC-2D pump (Eksigent™)
- Autosampler (Leap CTC™) equipped with 6-port micro-valve (VICI®) containing 1.0 μL loop
- Column #1 - PicoFrit® column (360 μm OD x 75 μm ID x 15 μm tip) packed with 20 cm ProteoPep II (5 μm, 300 Å C18, New Objective)
- Column #2 - PicoFrit column (360 μm OD x 75 μm ID x 15 μm tip) packed with 10 cm ProteoPep™ III (1.8 μm 90 Å C18, New Objective)
- Platinum microelectrode embedded in Teflon® elastomer (PicoClear™ conductive union, New Objective)

Reagents

- BSA digest (Waters®)
- 0.1% Formic Acid in water (JT Baker®)
- 0.1% Formic Acid in Acetonitrile (JT Baker)

Conditions

- Gradient: 8 or 20 minutes; 2%B to 50% B
- Flow rate: 300 or 500 nL/min.
- On-column injection: 300 fmol BSA digest

TABLE 1 - Changes in Column Pressure at Elevated Temperature

FlowRate (nL/min)	Temp (°C)	Column Pressure PSI (bar)	Temp (°C)	Column Pressure PSI (bar)	%Change
200	23	1750(121)	50	1350(93)	23%
300	23	2650(183)	50	1940(134)	27%
400	23	3446(238)	50	2530(174)	27%
500	23	4270(294)	50	3170(219)	26%

Column = 1.8 μm C18, 75 μm x 10 cm
Mobile phase = 2% ACN, 0.1% formic acid

FlowRate (nL/min)	Temp (°C)	Column Pressure PSI (bar)	Temp (°C)	Column Pressure PSI (bar)	%Change
200	23	686(47)	50	580(40)	15%
300	23	990(68)	50	838(58)	15%
400	23	1280(88)	50	1094(75)	15%
500	23	1575(109)	50	1365(94)	13%
600	23	1880(130)	50	1621(112)	14%
700	23	2169(150)	50	1884(130)	13%

Column = 3.5 μm C18, 75 μm x 10 cm
Mobile phase = 2% ACN, 0.1% formic acid

Recorded column pressure for two different columns at increasing flow rates. Column #1, packed with a 1.8 μm C18, exhibits an average 26% decrease in column pressure at an elevated temperature of 50°C as compared to 23°C. Column #2, packed with a 3.5 μm C18, exhibits an average 14% decrease in column pressure at an elevated temperature of 50°C.

